

CLOSED WOUND DRAIN SYSTEM (SPRING WITH PVC DRAIN)

- This product including 3-spring evacuator, PVC connecting tubing with Y connector, PVC round perforated drains with trocar.
- The design of 3-spring evacuator minimizes the Collapse-Resistance.

Item No.	Size	Item No.	Size
CWSP-9607	200ML+7FR	CWSP-9707	400ML+7FR
CWSP-9610	200ML+10FR	CWSP-9710	400ML+10FR
CWSP-9612	200ML+12FR	CWSP-9712	400ML+12FR
CWSP-9614	200ML+14FR	CWSP-9714	400ML+14FR
CWSP-9616	200ML+16FR	CWSP-9716	400ML+16FR
CWSP-9618	200ML+18FR	CWSP-9718	400ML+18FR



CLOSED WOUND DRAIN SYSTEM (SPRING WITH SILICONE DRAIN)

- This product including 3-spring evacuator, PVC connecting tubing with Y connector, Silicone round fluted drains with trocar.
- The design of 3-spring evacuator minimizes the Collapse-Resistance.

Item No.	Size	Item No.	Size
CWSS-9615	200ML+15FR	CWSS-9715	400ML+15FR
CWSS-9619	200ML+19FR	CWSS-9719	400ML+19FR

