

- This product includes silicone reservoir, silicone flat drains or other silicone drains and stainless steel trocar.
- Can use separately or using trocar for insertion and combining with silicone reservoir for fluids suction and collection.
- Provided with internal anti-reflux valve to prevent backflow of fluid to the patient.
- Silicone bulbs are suitable for plastic surgery, neurosurgery, orthopaedics, OB/GYN surgery and laparoscopic surgery

Item No.	Size ( ml )
----------	-------------

CWD-97XXXX	100,150,200,400
------------	-----------------

**Remarks:** Trocar , Drains and Reservoir can be Customised.

Item No.	Size
----------	------

Trocar	6inch , 8inch , Hook type , Bendable type , Double Bendable type.
--------	---

