



- Made of imported medical grade silicone.
- Flat four channels shape allows for greater multi-channel drainage.
- Channel drain offers important advantages over regular perforated drains by capillary effect.
- Radio-opaque line for easy X-ray visualization.
- Full fluted and 3/4 part fluted two types available.
- Available with or without “Three Face” stainless steel trocar.

Item No.	Size (mm)	Type
FF1-9707	7	3/4 Part Fluted
FF1-9710	10	3/4 Part Fluted
FF2-9707	7	Full Fluted
FF2-9710	10	Full Fluted