



- Made of imported medical grade silicone.
- Combines a transparent silicone round tube and white perforated flat drain.
- The inner wall of implanted flat segment with tooth-like structure facilitates drainage and lower the risk of occlusion.
- Radio-opaque line through the length for X-ray visualization.
- Full perforated and 3/4 part perforated two types available.
- Available with or without “Three Face” stainless steel trocar.

Item No.	Size ( mm )	Type
FP1-9707	7	3/4 Part Perforated
FP1-9710	10	3/4 Part Perforated
FP2-9707	7	Full Perforated
FP2-9710	10	Full Perforated