



- Made of imported medical grade silicone.
- Utilize internal channels or fluted, to swiftly remove exudate. Independent channels facilitate drainage and reduce risk of occlusion.
- Made by one time, hubless with no glue.
- Radio-opaque line through the length for X-ray visualization.
- Available with or without “Three Face” stainless steel trocar.

| Item No. | Size ( FR ) | with connector |
|----------|-------------|----------------|
| RF-9710  | 10FR        | XS             |
| RF-9715  | 15FR        | without        |
| RF-9719  | 19FR        | without        |
| RF-9724  | 24FR        | L              |
| RF-9728  | 28FR        | L              |