



- Made of imported medical grade silicone.
- Radio-opaque line through the length for X-ray visualization.
- Available with or without “Three Face” stainless steel trocar.

Item No.	Size ( FR )	with connector	Item No.	Size ( FR )	with connector
RP-9707	7FR	XS	RP-9716	16FR	without
RP-9710	10FR	XS	RP-9718	18FR	without
RP-9712	12FR	XS	RP-9720	20FR	without
RP-9714	14FR	S	RP-9722	22FR	L
RP-9715	15FR	without	RP-9724	24FR	L