

SILICONE RESERVOIR BOTTOM OUTLET

- Silicone bulbs are suitable for plastic surgery, neurosurgery, orthopaedics, OB/GYN surgery and laparoscopic surgery.
- Bottom outlet design are available for size 100ml and 200ml silicone reservoir.

Item No.	Size (mm)	Graduation (ml)
----------	-------------	-----------------

SRB-0910	100	25/50/75/100
SRB-0920	200	50/100/150/200



DRAINAGE BAG

- The drainage bag is designed to be connected with our silicone reservoir to collect, store and record the fluid after the surgeries.
- Be equipped with internal one-way valve to prevent reflux.
- Be equipped with an extension PVC tube and a luer lock to be connected with our silicone reservoir.
- Can be replaced individually.

Item No.	Size (ml)
----------	-------------

DB1000	1000
--------	------

